



TEST PROCEDURE

Check Instructions for Use for complete procedure



Result Interpretation:



CONTACT US

JOYSBIO (Tianjin) Biotechnology CO., LTD.


Headquarters: Tianjin International Joint Academy of
Biotechnology & Medicine 9th floor, No.220,
Dongting Road, TEDA 300457 Tianjin China

Website: <https://en.joysbio.com/>
Email: bd@joysbio.com
Tel: +86-186-0545-0986 (ENGLISH)
+86-22-5851-5062 (中文)
WhatsApp: +86-186-0545-0986

ORDER INFORMATION

JOYSBIO Monkeypox Virus Antigen Rapid Test Kit
(Colloidal Gold)
REF NUMBER:
G03821/G03822/G03823/G03824/G03825/G03826/G
03827/G03828
for 1/5/10/12/15/20/25/40 tests/box

© 2020 JOYSBIO (Tianjin) Biotechnology CO., LTD. All rights reserved.
All trademarks referenced are trademarks of JOYSBIO. Any photos
displayed are for illustrative purposes only.



Monkeypox Antigen Rapid Test Kit (Colloidal Gold)

Rapid monkeypox virus detection
with skin lesion exudate, acne
scab swabs or nasal swabs

JOYSBIO (Tianjin) Biotechnology CO., LTD.
Website: <https://en.joysbio.com/>



INTENDED USE

The Monkeypox Antigen Rapid Test Kit is a rapid, lateral flow immunoassay intended for the qualitative detection of monkeypox virus antigens from skin lesion exudate, acne scab swabs, or nasal swab specimens. This test is intended for use in individuals with or without symptoms or other epidemiological reasons to suspect a monkeypox infection. Self-test by a layperson is subject to local legislations.

HIGHLIGHT


- CE-IVD certified product
- Rapid detection of monkeypox virus antigen in less than 15 minutes
- No instrumentation required
- Cost-effective solution for large-scale screening

SPECIFICATION:

- **TEST TIME:** 15-20 MINUTES
- **STORAGE CONDITION:** 2°C-30°C
- **CE-IVD MARK**
- **SAMPLE TYPE:** Swabs
- **SHELF LIFE:** 24 months

CONTENTS:

- Test Device
- Buffer Bottle (350 µL/bottle)
- Extraction Tube with Nozzle Cap
- Sterilized Swab
- Instructions for Use



PACKING INFO

Package	# of boxes/test s per carton	Gross weight /volumetric weight (kg)
1 tests/box	472/472	7.5/21.5
5 tests/box	219/1095	17.39/21.5
20 tests/box	50/1000	15.88/21.5
25 tests/box	50/1250	18.39/21.5

Carton size: 410x510x610mm
Carton volumetric weight: 21.5kg

